

# 33 Gal Red Biohazard Waste Bag, Pack of 100 - 1691

## Product Images







## Short Description

---

### Key Benefits

- Specifically labeled for Biohazard material handling
- Linear Low Density Polyethylene (LLDPE) for improved flexibility and shear resistance
- Fits Eagle 30 Gallon Nestable Drums

## Description

---

For holding or consolidating larger amounts of biological waste or contaminated materials. 33 Gallon capacity bag is designed to fit Eagle 1601 series drums. Once the bag is full and ready for disposal, simply seal the bag in place and close the drum using its lid and lever-lock. Red color bag is easily identifiable as biohazard bags, drums can be Black, Yellow, Blue or Red. Printing includes English, Spanish and French translation along with the international recognizable biohazard symbol. Using Linear Low Density Polyethylene (LLDPE) material at 1.2 Mil thickness, a high level of flexibility and shear resistance can be attained without leaks. Meets ASTM D 1709-01 and ASTM D 1922-00a specifications. A tubular construction, single side seal and star seal bottom provide a stable, leak proof bag. Sold in packs of 100.

## Specifications

---

Model No	1691
Color	Red
UPC	048441693023
Gallon Capacity	33
Liter Capacity	124
Material Specifications	Linear Low-density Polyethylene
Net Weight, lbs	1.17
Net Weight, kgs	0.53
Wall Thickness	1.2 Mil
Single Bag Size	32.5" H x 32.5" W x 45" D