

Under Fume Hood Flammable Cabinets

Product Images





Short Description

- **Durable Construction:** Built with welded, 18-gauge steel double-walls and 1.5 inches of insulating air space for superior fire resistance
- **Versatile Options:** Available in capacities from 15 to 31 gallons and various color options
- **Secure Closure:** Features a stainless steel, 3-point bullet self-latching system for secure closure, available in manual closing or self-closing options
- **Spill Management:** Includes adjustable SpillSlope shelves directing spills to a leakproof 2-inch sump, ensuring safety and compliance
- **Compliance and Warranty:** Meets OSHA and NFPA safety regulations and includes a ten-year warranty; FM Approved
- **Ergonomic Design:** Features a 4-inch high by 3-inch deep toe kick for safe working distance and measures 35.75 inches in height to fit under fume hoods or counters

Description

Justrite Under Fume Hood Safety Cabinets are designed to safely store solvents and flammable liquids in laboratory settings, ensuring compliance with safety regulations and enhancing workplace safety. These cabinets are constructed with welded, 18-gauge steel double-walls and 1.5 inches of insulating air space, providing superior fire resistance. They feature back wall vents with built-in flame arresters and a stainless steel, 3-point bullet self-latching system for secure closure.

All under fume hood safety cabinets measure 35.75 inches in height, making them ideal for fitting under fume hoods, benches, or counters in laboratories. The ergonomic design includes a 4-inch high by 3-inch deep toe kick, allowing users to maintain a safe distance while working. The cabinets feature adjustable SpillSlope shelves that direct spills to a leakproof 2-inch sump, ensuring safety and compliance. A durable powder-coat finish provides maximum chemical resistance inside and out.

To configure your under fume hood safety cabinet, select from various color options, including yellow, silver, or White/light neutral, to complement your laboratory setting. Choose from capacities ranging from 15 to 31 gallons to meet your specific storage needs. For the 15-gallon capacity, you will also have the option to select the door hinge side. Additionally, choose between manual closing or self-closing doors to meet your specific safety and compliance requirements.

Self-closing doors are a critical safety feature in these cabinets, offering enhanced protection against fires and accidental releases. These doors automatically close and latch in the event of a fire or emergency, using fusible links that melt at 165°F (74°C), causing the doors to shut automatically. This feature minimizes human error by ensuring that the cabinet is always properly closed, even when staff are distracted or forgetful. While manual closing doors are available, self-closing doors provide the safest option, especially in laboratory environments where safety is paramount.

Justrite Under Fume Hood Safety Cabinets are designed to meet the complex requirements of handling and storing hazardous materials in laboratories, providing a reliable solution for safe storage and compliance. By configuring your cabinet, you ensure a safe and compliant storage solution tailored to your laboratory's unique needs.

All Justrite cabinets carry a ten-year warranty and are FM-approved, ensuring peace of mind and reliable performance.

Specifications

Model No	88-UH
----------	-------

Prop 65

⚠️WARNING: This product can expose you to chemicals including lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Product Options

Color:	Yellow
	Light Neutral
	Silver
Door Style:	Self Close
	Manual
Capacity:	15 gallon
	19 gallon
	23 gallon
	31 gallon
Cabinet Door Hinge:	Right
	Left
	None